




THE LUTING COMPOSITE FOR
ALL INDIRECT
RESTORATIONS



ON THE SAFE SIDE –
VISALYS® CEMCORE

**DUAL-CURING LUTING COMPOSITE
AND CORE BUILD-UP MATERIAL**



KETTENBACHDENTAL
Simply intelligent

VISALYS® CEMCORE – ENTIRELY CONVINCING

FOR ALL CASES

- **Suitable for permanent cementation** of crowns, bridges, inlays, onlays, veneers and adhesive bridges, root posts and indirect core build-ups
- **Suitable for all dental restoration materials** ((noble)metal, oxide ceramic, (zirconium oxide and aluminum oxide ceramics), silicate ceramic (feldspar and glass ceramic), hybrid ceramic, composite)

Conclusion by Prof. Dr. Claus-Peter Ernst

[...] (One) impressive potential by **Visalys® CemCore** for a huge range of indications [...]

Source: Ernst, Prof. Dr. Claus-Peter et al.: A universal adhesive cementation and core build-up material, in: Spitta: ZMK, Vol. 36, Issue 5 May 2020 [294–308]



HIGH COLOR STABILITY

Unlike other luting composites, **Visalys® CemCore** is hydrophobic. This minimizes the absorption of colors and discoloration of the cement joint.



	RelyX™ Universal	RelyX™ Unicem	RelyX™ Ultimate	Panavia™ F2.0	Panavia™ V5	Multilink® Automix	Variolink® Esthetic DC	Visalys® CemCore
Water								
Coffee								
Red wine								

Light-cured samples, stored 48 h in the corresponding solution at 37 °C.

SECURE ADHESION & LONG LIFE




The **high quality of the marginal gap** of Visalys® CemCore represents the high quality of the bond between the dental substance and the composite.

Source: Prof. Dr. Roland Frankenger et al. (2019), Marginal Quality of Ceramic Inlays luted with Visalys® CemCore, University Marburg

No debonding and no fracture, even after **10 years** of simulated oral use.

Source: Prof. Dr. Dipl.-Ing. Martin Rosentritt (2020), Fracture force after long-term storage and TCML, University Regensburg



SMART HANDLING

Always the same procedure, regardless of the type of indirect restoration.



Easy to remove any excess

1) Outstanding adhesion results with patented Active-Connect-Technology.

More Information: www.kettenbach-dental.com/act

2) Copyright: Dr. med. dent. Frank-Michael Maier



VISALYS® CEMCORE – THE RANGE

Visalys® CemCore

Starter pack	4.5 g (2.5 mL) automix syringe, 2 mL Tooth Primer, 2 mL Restorative Primer, blue mixing tips, blunt/tapered, 5 units each, 3 intraoral tips, 2 endo tips	Universal (A2/A3)	REF 13570
Normal pack	9 g (5 mL) automix syringe, blue mixing tips, blunt/tapered, 10 units each, 6 intraoral tips, 4 endo tips	Universal (A2/A3)	REF 13572
		Opaque	REF 13573
		Translucent	REF 13574
		Bleach	REF 13575
		Dark (A4)	REF 13576
Visalys® Tooth Primer	4 mL bottle		REF 13580
Visalys® Restorative Primer	4 mL bottle		REF 13581



Made
in
Germany



1800 651 788 www.pharmacodental.com.au

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE.

Some medical devices are for healthcare professional use only.
Some products or variants may not be registered or available in Australia.