

PREVENTION

RESTORATION

IMPRESSIONS

SILGINAT® – FOR EVEN GREATER EFFICIENCY!

A-SILICONE FOR ANATOMIC IMPRESSIONS

UNIQUE: 100% DEVELOPED AND MANUFACTURED IN GERMANY

We oversee every step – from initial concept to final product!
True German manufacturing, true responsibility, and quality you can truly trust.



MADE IN
GERMANY




POURABLE MULTIPLE TIMES,
ONLY ONE ANATOMICAL
IMPRESSION PER PATIENT
REQUIRED.



KETTENBACHDENTAL
Simply intelligent

TIME SAVINGS AND REDUCED WORKLOAD

INDICATIONS



Technical data

	Working time at 23°C ≤ / 74 °F ≤	Intraoral setting time at 35°C ≥ / 95 °F ≥	Total setting time* ≥	Elastic recovery test (approx.) %	Linear dimensional change (approx.) %	Hardness (approx.) Shore
 Silginat® / Silginat® Strawberry 1:1 cartridge	1:15 min.	1:15 min.	2:30 min.	99.5	– 0.20	A 42
 Silginat® / Silginat® Strawberry 5:1 cartridge	1:30 min.	1:30 min.	3:00 min.	99.5	– 0.20	A 45

*Total setting time (removal time from mouth) from the start of mixing

SILGINAT® – READY FOR PRACTICAL TESTING

Simply more options. Find your ideal product now!



Anatomic impressions

SILGINAT® Normal pack	6 x 50 mL cartridges, 6 mixing tips, white	REF 13846
SILGINAT® Bulk pack	24 x 50 mL cartridges	REF 13847
SILGINAT® Refill pack	2 x 380 mL cartridges	REF 14713
SILGINAT® Strawberry Normal pack	6 x 50 mL cartridges, 6 mixing tips, white	REF 13826
SILGINAT® Strawberry Bulk pack	24 x 50 mL cartridges	REF 13827
SILGINAT® Strawberry Refill pack	2 x 380 mL cartridges	REF 14715



Accessories

Universal Adhesive	10 mL bottle	REF 14206
Dynamic mixers	45 units, white	REF 17900
Mixing tips	50 units, white 100 units, white	REF 17226 REF 17240
Sympress Dispenser	Electrical dosing and mixing devices	REF 35910
Applyfix 4	For 50 mL cartridges 1:1	REF 17204



1800 651 788 www.pharmacodental.com.au

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE.

Some medical devices are for healthcare professional use only.
Some products or variants may not be registered or available in Australia.