

**STRONGER, GLOSSIER AND MORE RELIABLE –
A NEW GENERATION OF PROVEN QUALITY!**



TEMPORARY CROWN AND BRIDGE MATERIAL



THE NEW GENERATION OF VISALYS® TEMP

VISALYS® TEMP ONE:ONE offers maximum stability, high fracture strength, and excellent edge stability. As a 1:1 crown and bridge material, it delivers precise, reliable, and aesthetically high-quality results – without compromise.



Minimal waste, maximum efficiency. Ideal for modern practices with the highest standards.

Bleed cartridge only once
**MULTIPLE APPLICATIONS
POSSIBLE***



Up to **62 ANTERIOR CROWNS**

SAVE 35% MORE on materials

**High
fracture strength
& edge stability.**

**Low
smear layer:**
Easier handling
after removal.

High gloss
without polishing.

**1:1
SYSTEM**

Precise milling:
Low dust generation,
precise results.


Easy to repair:
Easy to supplement
after curing.




RESTAURATION

... from **TEMPORARY** to **PERMANENT** – uncompromising strength

OUTSTANDING FRACTURE STRENGTH



CONVENIENT WORKING AND SETTING TIMES



Note: The working times apply to products stored and processed at a temperature of 23°C/74°F and 50 % relative humidity. Higher temperatures reduce the working time and lower temperatures extend these times.

VISALYS® TEMP ONE:ONE
GET TO KNOW

VISALYS® TEMP ONE:ONE – PRODUCT RANGE

Color	Normal pack	Bonus pack
	1 x 50 mL (73 g) cartridges, 15 orange mixing tips	5 x 50 mL (73 g) cartridges, 15 orange mixing tips
Bright (B1-A1)	REF 13850	
Medium (A2-A3)	REF 13851	REF 13856
Natural Bleach	REF 13849	
Bleach	REF 13853	
Mixing tips orange, 60 units		REF 17217



EXTRAS



DREAMTEAM

Silginat® Refill pack	2 × 380 mL cartridges	REF 14713
Dynamic mixers	45 units	REF 17900
Multi Tray Anterior	30 units	REF 17752
Multi Tray Posterior	50 units	REF 17750
Multi Tray Quadrant	30 units	REF 17753
Orthoskavident® C Normal pack	1 × 150 mL bottles	REF 13063
Orthoskavident® C Economy pack	3 × 150 mL bottles	REF 13065
Applyfix 4 (1:1)		REF 17203